


The following can be printed and used as reference and for recording foot measurements, which are then entered into LVP4.0 using the FOOTC data entry page. Measurements should be entered in mm. Measurements are most accurately taken using a digital caliper but any suitable method may be used. If you are interested in obtaining a suitable digital caliper contact us.

|                |       |
|----------------|-------|
| Ly             | _____ |
| Lxy            | _____ |
| W <sub>x</sub> | _____ |
| W <sub>z</sub> | _____ |
| W <sub>y</sub> | _____ |
| L12            | _____ |
| L8             | _____ |
| W12            | _____ |
| W8             | _____ |
| W4             | _____ |
| H <sub>x</sub> | _____ |
| H <sub>z</sub> | _____ |
| H <sub>y</sub> | _____ |



For ease of data entry and consistency, enter data into the FOOTC page of LVP4.0 in millimeters (mm) only

**This page can be printed and used as a Data Recording Sheet for Your foot measurements**

Copyright 2006 BIMECO-BSRI [support@limbvolumes.org](mailto:support@limbvolumes.org)